

Pierre Aller, Matthew Rodrigues and Abbey Telfer

XFEL-Hub @ Diamond light source

pierre.aller@diamond.ac.uk

matthew.rodrigues@diamond.ac.uk

a.telfer@diamond.ac.uk


XFEL-Hub & time-resolved serial crystallography

~~17th September 2025~~




Overview:

- XFEL-Hub introduction
- X-ray sources: Synchrotron and XFEL
- Sample deliveries current and upcoming
 - High viscosity extruder
 - Fixed target
 - Hybrid method: drop on tape

- 
- XFEL-Hub introduction
 - X-ray sources: Synchrotron and XFEL
 - Sample deliveries current and upcoming
 - High viscosity extruder
 - Fixed target
 - Hybrid method: drop on tape



XFEL Hub



Allen Orville, Principal Scientist XFEL-Hub

Pierre Aller, Senior Beamline Scientist

Matthew Rodrigues, Beamline Scientist

Abbey Telfer, PDRA

Tiankun Zhou, PDRA

Emily Freeman, PhD student

Jack Stubbs, PhD student

Eloisa Wheatley, PhD student

Madeleine Smith, year in industry student

<http://www.diamond.ac.uk/Beamlines/Mx/XFEL-Hub.html>

In This Section

XFEL Travel Assistance

Deadlines

Staff

Sample delivery methods

Dynamic Structural Biology BAG at Diamond

Gallery

XFEL-Hub contact

XFEL-Hub@diamond.ac.uk

XFEL-Hub Principal Scientist:

Allen Orville

Tel: +44 (0) 1235 567505

E-mail: allen.orville@diamond.ac.uk

Welcome to the XFEL Hub at Diamond


The UK is taking a leading role in the development of a new structural biology facility (SFX) at the European X-ray Free Electron Laser (XFEL), in Hamburg, Germany, and a complementary facility at Diamond (The UK XFEL Hub) to help develop the required expertise.

The UK is a leading member of the serial femtosecond crystallography user consortium (SFX-UC), and through its role in the collaboration will deliver a dual instrument at the European XFEL that will enable visiting scientific users to image macromolecules from nanocrystals on the SFX instrument and single particles using scattering of coherent X-rays and detecting the coherent diffraction pattern on the [SPB instrument](#).



New MFX end station at LCLS, SLAC National Accelerator Laboratory, USA

Instruments by Science Group



Serial femtosecond crystallography (SFX) was the term coined to describe the protocols used to collect the first diffraction data from nanocrystals at the [LCLS in Stanford, USA](#). Scientists there were able to exploit the ultrashort pulses generated by LCLS to collect single shot diffraction images from nanocrystals. Diffraction from such small crystals was made possible due to the high brilliance of the beam (many orders of magnitude more than achievable at a synchrotron).

The SFX-UC brings together members of the scientific communities in Germany, United Kingdom, Sweden, Slovakia, Switzerland, United States and Australia. The project is being managed in Hamburg by Adrian Mancuso. The UK structural biology community is represented by Jim Naismith (lead PI, St. Andrews), Jan Lowe (LMB, Cambridge), Henry Chapman (CFEL, Hamburg), Martin Walsh (Diamond) and Dave Stuart (Oxford & Diamond).

At Diamond Dr Allen M. Orville has been appointed group leader of the XFEL Hub which will act as a focus for a number of activities:

- Development of hardware and software for SFX
- Provision of a sample environments lab to enable users to prepare samples for current liquid jet technologies at operating hard X-ray FELS
- Provision of a user access program for SFX/SPB and currently operating hard X-ray FEL facilities worldwide.

Housed within the existing Diamond infrastructure, the hub will enable users to fully prepare for their experiments with currently operating XFELs and the European XFEL when it comes online in Hamburg in 2017. The UK Hub will provide support in terms of sample preparation, data processing and training. There will also be a dedicated fibre link from Hamburg to Harwell enabling users to carry out data analysis back in the UK, with support from the UK Hub team.

For more information on the UK-XFEL Hub, please contact [Allen Orville](mailto:Allen.Orville).

More about SFX and the European XFEL see [latest news](#).

XFEL-Hub at Diamond

Travel grant for XFEL experiment:


- Only UK based researchers
- Up to 4 members on the same experiment can apply
- Non-proprietary XFEL beamtime
- Only flight ticket will be reimbursed



Dynamic Structural Biology BAG

- Open to anyone interested in SSX
- ~ 10 shifts available on I24 /Allocation period
- Sample delivery systems available: Fixed target and LCP extruder



- 
- XFEL-Hub introduction
 - X-ray sources: Synchrotron and XFEL**
 - Sample deliveries current and upcoming
 - High viscosity extruder
 - Fixed target
 - Hybrid method: drop on tape

Synchrotron facilities



emittance: measure of size and angular spread of the e^- beam

- 1st Generation sources
 - parasitic
- 2nd Generation sources (80's)
 - dedicated machines to synchrotron rad. emittance typically high (150nmrad)
- 3rd Generation sources (90's-00's)
 - low beam emittance, use of insertion devices (wiggler and undulators) – emittance on order of 3-20nmrad
- 4th Generation sources (20's)
 - Change the lattice, lower beam emittance ~ 130 pmrad



Key / Légende

	Operational 3rd generation synchrotron Synchrotron de 3ème génération opérationnel
	3rd generation synchrotron under construction Synchrotron de 3ème génération en construction
	Operational 2nd generation synchrotron Synchrotron de 2ème génération opérationnel
	Proposed new synchrotron Nouveau synchrotron proposé

Image from Canadian Light Source

XFEL sources

Linear accelerator followed by undulators
- 700-1.1km (Sacla, PalXFEL, SwissFEL)
- 3km LCLS-II
- 3.4km European XFEL



SACLA



European XFEL



XFEL facilities



5 XFEL facilities:

- **LCLS-II** (USA)
- **SACLA** (Japan)
- **PAL XFEL**
(Korea)
- **SwissFEL**
(Switzerland)
- **European XFEL**
(Germany)




Synchrotron vs XFEL: Beam properties

XFEL	Synchrotron
fs pulse	Continuous (ish)
Coherent beam	Incoherent beam
Limited access	Easy access

$\sim 10^{10}$ more brilliant

Photons in 1 XFEL pulse \simeq 1s of synchrotron beam



From Shintake T., Proceedings of PAC07,
Albuquerque, New Mexico, USA




Synchrotron vs XFEL



Slow shutter speed

Synchrotron




Faster shutter speed


XFEL

Shutter always open
Very short flash of light

=> Better temporal definition



Motionless object



Dynamic structural biology



Intermediate states



Resting state at
 $t=0$




Final state



Falling Cat, 1894 by Etienne-Jules Marey

diamond

Synchrotron vs XFEL: Time domain resolution



Orville, A. M. (2020)
Current Opinion in Structural Biology
65, 193-208

XFEL

X-ray Free Electron Lasers

Diamond II

Storage ring (synchrotron)

Diamond

Diffraction-limited storage ring


MeV ED

cryo ED

Electron diffraction (RT / cryo)

cryo EM

Single particle cryo-electron microscopy




XFEL


Endstation

Synchrotron

User friendly




MFX experimental hutch LCLS



I24 Diamond




- 
- XFEL-Hub introduction
 - X-ray sources: Synchrotron and XFEL
 - Sample deliveries current and upcoming**
 - High viscosity extruder
 - Fixed target
 - Hybrid method: drop on tape


Sample delivery methods available for SSX

Jet

GDVN, Co-Mesh, mixing injector, LCP extruder...



Owen et al., Acta D 2017



Fixed target



I24's portable station,
Road runner, goniometer...



Hybrid method

Acoustic droplet ejection (ADE) / drop
on tape

High Viscosity Extruder




- Low sample consumption <0.3 $\mu\text{l}/\text{min}$
- Crystal size up to 30 – 50 μm
- Can be used at both Synchrotron and XFEL
- Designed for membrane protein in LCP originally
- Install on I24 upon request
- Pump probe experiment

Weierstall et al., Nat Commun (2014)


Zhou et al., Adv Exp Med Biol (2016)

Fixed target: with a chip

Crystal(s) fixed on sample holder and during data collection the sample holder move to different positions



- <100 µl Crystal solution pipetted over chip made of silicon nitride
- Chip capacity = 25,600
- Funnel-like apertures for holding crystals



- low consumption
- Up to 40 µm crystals
- For XFEL and synchrotron
- Time resolved pump probe and mixing experiment on i24



Drop on Chip (DoC) mixing experiment for time resolved serial crystallography

nature methods

SUPPLEMENTARY INFORMATION

<https://doi.org/10.1038/s41592-019-0553-1>

In the format provided by the authors and unedited.

Liquid application method for time-resolved analyses by serial synchrotron crystallography

Pedram Mehrabi^{1,7}, Eike C. Schulz^{1,2,7}, Michael Agthe³, Sam Horrell³, Gleb Bourenkov⁴, David von Stetten⁴, Jan-Philipp Leimkoh⁵, Hendrik Schikora⁵, Thomas R. Schneider⁴, Arwen R. Pearson^{2,3}, Friedjof Tellkamp⁵ and R. J. Dwayne Miller^{1,2,6*}

¹Max-Planck-Institute for Structure and Dynamics of Matter, Department for Atomically Resolved Dynamics, Hamburg, Germany. ²Centre for Ultrafast Imaging, Universität Hamburg, Hamburg, Germany. ³Institute for Nanostructure and Solid State Physics, Universität Hamburg, Hamburg, Germany.


⁴European Molecular Biology Laboratory, Hamburg Unit c/o Deutsches Elektronen Synchrotron, Hamburg, Germany. ⁵Max-Planck-Institute for Structure and Dynamics of Matter, Scientific Support Unit Machine Physics, Hamburg, Germany. ⁶University of Toronto, Departments of Chemistry and Physics, Toronto, Ontario, Canada. ⁷These authors contributed equally: Pedram Mehrabi, Eike C. Schulz. *e-mail: dwayne.miller@mpsd.mpg.de




Autodrop pipette
70 µm aperture
~ 180-200 pL
Up to 6kHz



a





b



Set up used at the T-REXX
beamline (Hamburg)

diamond

The logo for Diamond Light Source, featuring a stylized yellow starburst icon followed by the word "diamond" in a lowercase sans-serif font.



Aperture cross-section
(x2000 magnification)

Apo state

Microcrystal


Substrate ejection

Substrate


Diffraction Δt

X-rays

I24 set-up



DoC at Diamond i24




Ejection of substrate on chip


>95 % hit

30 µm ID capillary ~ 50-70 pl drop


10 ms dwell time at each well

About 75 time slower than real time





'+ ligand' dataset
 $F_o - F_o^{Apo} (\pm 3\sigma)$



Lysozyme \pm GlcNAc, $\Delta t = 1900$ ms

control dataset
 $F_o - F_o^{Apo} (\pm 3\sigma)$



1 chip 2 datasets


- Interleaved control
- Substrate added

Possible leakage between apertures
 \Rightarrow unreliable Δt

Sample delivery methods available for SSX

Jet


GDVN, Co-Mesh, mixing injector, LCP extruder...



Owen et al., Acta D 2017

Fixed target

I24's portable station,
Road runner, goniometer...




Hybrid method

Acoustic droplet ejection (ADE)



On demand: drop on tape

hybrid delivery method using both jet (acoustic droplet ejection) and fixed target (tape drive)




Some Characteristics of Acoustic Droplet Ejection (ADE)

Transducer frequency (MHz)	~ Droplet volume (nL)	~ Droplet diameter (μm)	~ Crystal size within drop (μm^3)
11.5	2.5	170	10 – 100
15	1	125	10 – 75
25	0.1	60	10 – 40

Max ejection frequency ~ 100 Hz

The ADE-tape drive sample delivery system for tr-SFX ± tr-XES



Fuller and Gul *et al.*, (2017) Nature Methods 14: 443-449


Butrym, Simon, Aller *et al.*, (2021) Nat Commun. 12(1):4461



Reaction initiation in crystals:



- a) Pump probe: laser illumination
- b) Mixing experiment:
 O_2 mixing in a reaction chamber
- c) Drop on drop induces turbulent mixing

- Various crystal size up to 100 μm
- Wide range of experiments
(modularity)
- Low consumption few $\mu l/min$
- For XFEL and synchrotron



Drop on drop induces turbulent mixing

In 2019 first test of mixing experiment at SACLA



Lysozyme + GlcNAc



Mixing experiment at SACLA and LCLS, $\Delta t = 10 \text{ ms} - 2\text{s}$

From Butryn, Simon, Aller *et al.*, Nat Commun. 2021





Subscribed


ARTICLE | November 8, 2023

In Situ Structural Observation of a Substrate- and Peroxide-Bound High-Spin Ferric-Hydroperoxo Intermediate in the P450 Enzyme CYP121

Romie C. Nguyen, Ian Davis, Medhanjali Dasgupta, Yifan Wang, Philipp S. Simon, Agata Butryn, Hiroki Makita, Isabel Bogacz, Kednerlin Dornevil, Pierre Aller, Asmit Bhowmick, Ruchira Chatterjee, In-Sik Kim, Tiankun Zhou, Derek Mendez, Daniel W. Paley, Franklin Fuller, Roberto Alonso Mori, Alexander Batyuk, Nicholas K. Sauter, Aaron S. Brewster, Allen M. Orville, Vittal K. Yachandra, Junko Yano, Jan F. Kern*, and Aimin Liu*




X-ray emission spectroscopy



Applications:



- Probing metal oxidation state for transition metals ($K\alpha_{1,2}$)
- Metal oxidation and spin state ($K\beta'_{1,3}$)
- Information about ligands and coordination sphere ($K\beta''_{2,5}$)



Adapted from
DeBeer & Bergmann (2016)
Encyclopedia of Inorganic and Bioinorganic Chemistry



Von Hamos Spectrometer



Current scope:

Fe	²⁶
iron	
$K\alpha_{1,2}$ 6405, 6392 eV	
$K\beta_{1,3}$ 7059, 7059 eV	
$K\beta_{2,5}$ 7110 eV	

Cu	²⁹
copper	
$K\alpha_{1,2}$ 8046, 8026 eV	
$K\beta_{1,3}$ 8903, 8901 eV	
$K\beta_{2,5}$ 8974 eV	

Ni	²⁸
<td></td>	
$K\alpha_{1,2}$ 7480, 7463 eV	
$K\beta_{1,3}$ 8266, 8265 eV	
$K\beta_{2,5}$ 8329 eV	

Future scope:


Mn	²⁵
manganese	
$K\alpha_{1,2}$ 5900, 5889 eV	
$K\beta_{1,3}$ 6491 eV	
$K\beta_{2,5}$ 6537 eV	

Co	²⁷
cobalt	
$K\alpha_{1,2}$ 6930, 6915 eV	
$K\beta_{1,3}$ 7649, 7650 eV	
$K\beta_{2,5}$ 7706 eV	


Mo	⁴²
molybdenum	
$K\alpha_{1,2}$ 17480, 17375 eV	
$K\beta_{1,3}$ 19606, 19588 eV	
$K\beta_{2,5}$ 19962, 19768 eV	

XES + XRD data collected from AcNIR crystals to detect oxidation state of Cu

Copper nitrite reductases are key enzymes in the nitrogen cycle.




Apo-state
 $R_{\text{free}} = 0.26$
 $R_{\text{work}} = 0.22$




The denitrification pathway in these organisms couples ATP synthesis with the reduction of nitrite (NO_2^-) to dinitrogen (N_2), via nitric oxide (NO) intermediates. This pathway is critical in the return of fixed nitrogen to the atmosphere and in controlling the level of biologically available nitrogen in the soil.

Dataset name	xia2-dials Integrated	xia2-dials- ransac Integrated	xia2-dials Merged	xia2-dials- ransac Merged	xia2-dials $\text{CC1/2} = 0.3$	xia2-dials- ransac $\text{CC1/2} = 0.3$
Apo-state 0.58 Mgy (reduced)	17755	20947	17748	20893	1.9	1.9
NO_2 bound 0.58 Mgy (reduced)	10920	14028	10893	13873	1.94	1.92


Recent Structural Insights into the Function of Copper Nitrite Reductases (2017)




Acnir + LCP using HVE injector



Dataset 412-415



Dataset 420-424





Acnir (Cu I)



Acnir (Cu II)





In conclusion:

⇒ How to choose your sample delivery?

Do I need an XFEL or a synchrotron source?

What do I want to achieve?

How much sample can I produce?

What is the crystal size and density?

If you want to try SSX, please contact us! We can offer beamtime at Diamond through our BAG!



Acknowledgements

Diamond Light Source

XFEL-Hub team

Allen Orville

Matt Rodrigues

Tiankun Zhou

Abbey Telfer

Eloisa Wheatley

Jack Stubbs

PolyPico

Gabriel Leen

University of Bristol

Jim Spencer

Phil Hinchliffe

Laura Parkinson

University of Bath

Cat Tooke

Joe Hoff

University of Leicester

Hanna Kwon

Abi Balakrishnan

Eu XFEL staff

Financial support, in part, to AMO:

- Wellcome Investigator Award in Science
(210734/Z/18/Z)

